**Multi-parameter Analyzer
a-AB33M1
Bench pH Meter
a-AB33PH
Bench Conductivity Meter
a-AB33EC
Instruction Manual**

Change records

Date	Version	Description
2024/9/10	A	New release

Contents

1.	Introduction	3
1.1.	Safety precautions	3
1.2.	Intended use	4
2.	Installation and setup	5
2.1.	Installation of stand-alone electrode holder	5
2.2.	Installation the power adapter	6
2.3.	Connect the pH and conductivity electrodes	6
2.4.	Touch screen	7
2.5.	Buttons	9
2.6.	Turn on and off the meter/analyzer	10
3.	User management	11
3.1.	User management	11
4.	Records management	12
5.	Setup menu	13
5.1.	Setup	13
5.1.1	Mode	14
5.1.2	pH (a-AB33M1, a-AB33PH)	16
5.1.3	Conductivity channel setting (a-AB33M1, a-AB33EC)	17
5.1.4	Temperature	18
5.2.	i-Steward™ menu	19
5.2.1	Calibration due alarm (Cal. due alarm)	19
5.2.2	pH ECS	19
5.2.3	Export log	20
5.2.4	Log clear (ID0 only)	21
5.2.5	Self-test	21
5.2.6	Reset	21
5.3.	Settings menu	22
5.3.1	Language	23
5.3.2	RS232	23
5.3.3	Date/time	24
5.3.4	Brightness	25
5.3.5	Keypad beep	25
5.3.6	Auto shut off	25
5.3.7	About meter	26
6.	pH operation (a-AB33M1, a-AB33PH)	27
6.1.	pH calibration process	27
6.1.1	pH electrode direct calibration process	27

6.1.2	ORP electrode (relative mV mode) calibration process.....	28
6.2.	pH measurement process	29
6.2.1	pH measurement.....	29
6.2.2	mV and ORP measurement	30
6.2.3	Temperature measurement.....	30
7.	Conductivity operation (a-AB33M1 and a-AB33EC).....	31
7.1.	Conductivity calibration process	31
7.2.	Conductivity measurement.....	31
7.3.	TDS, SAL (Salinity) and RES (resistivity) measurement.....	32
8.	Print	33
8.1.	Format	33
8.2.	Output format.....	36
9.	Commands	36
10.	Maintenance	37
10.1.	Error message.....	37
10.2.	Meter maintenance	38
10.3.	Electrode maintenance	38
10.3.1	pH electrode maintenance.....	38
10.3.2	Conductivity electrode maintenance	38
10.4.	Cleaning	38
10.5.	Customer service.....	39
11.	Technical data.....	40
11.1.	Specification	40
11.2.	Size	45
12.	Appendixes.....	46
12.1.	pH buffer group	46
12.2.	Conductivity standards.....	47
12.3.	Example of Temp. Coeff. (α value)	48
13.	Compliance.....	49

1. Introduction

This manual contains installation, operation and maintenance instructions of a-AB33M1 bench multi-parameter analyzer, a-AB33PH bench acidimeter and a-AB33EC bench conductivity meter. Please read this manual carefully before installation and operation.

The table below lists the measurement modes each model supports:

Table 1-1 Measurement modes

Model	Supported measurement modes
a-AB33PH	pH, mV, relative mV and ORP measurement
a-AB33EC	Conductivity, TDS, salinity and resistivity measurement
a-AB33M1	pH, mV, relative mV and ORP measurement on channel one; Conductivity, TDS, salinity and resistivity measurement on channel two

1.1. Safety precautions

Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the meter/analyzer, malfunctions and false results.

- WARNING** For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
- CAUTION** For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
- ATTENTION** For important information about the product. May lead to equipment damage if not avoided.
- NOTE** For useful information about the product.

Warning Symbols

General hazard

Explosion hazard

Electrical shock hazard

Safety precautions

CAUTION: Read all safety warnings before installing, making connections or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

- Before connecting power, verify that the AC adapter's input voltage range and plug type are compatible with the local AC mains power supply.
- Do not position the equipment such that it is difficult to reach the power connection.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- The equipment is for indoor use only.
- Do not operate the equipment in wet, hazardous or unstable environments.
- When using chemicals and solvents, comply with the instructions of the chemical producer and the general lab safety rules.
- Do not allow liquids to enter the equipment.
- Disconnect the equipment from the power supply when cleaning.
- Service should only be performed by authorized personnel.

WARNING: When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules.

WARNING: Never work in an environment subject to explosion hazards! The housing of the meter/analyzer is not gas tight. (Explosion hazard due to spark formation, corrosion caused by the ingress of gases).

WARNING: Electrical shock hazards exist within the housing. The housing should only be opened by authorized and qualified personnel. Remove all power connections to the unit before opening.


1.2. Intended use

This meter/analyzer is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for measuring the parameters described in these instruction manuals. Any other type of use and operation beyond the limits of technical specifications, without written consent from OHAUS, is considered as not intended. This meter/analyzer complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the meter/analyzer is not used according to these operating instructions, the intended protection provided by the meter/analyzer may be impaired.


2. Installation and setup

2.1. Installation of stand-alone electrode holder

1. Attach the magnetic part on the bottom of the meter to the magnetic part on the stand-alone electrode holder to finish installation. (Both marked in red demonstrated below.)

2. The two tension knobs can be adjusted by turning clockwise to tighten or turning counterclockwise to loosen. When they are loosened, you can move the holder up and down, forward and backward to the desired angle. Then you can tighten the knobs to fix the holder.

2.2. Installation the power adapter

The Meter/Analyzer comes with a universal power adapter, so as to ensure that an appropriate power adapter is connected to your power supply. Using other power adapters will damage the meter/analyzer and void the warranty.

- 1 Plug the correct connector plug into the AC adapter until it is fully plugged.
- 2 Connect the AC adapter wire and the DC slot of this meter/analyzer.
- 3 Install the cable in this manner to ensure that it will not be damaged or interfere with the operation.
- 4 Plug the plug of the AC adapter into an accessible power outlet. In addition, it is recommended to use a surge protection device (SPD) or uninterrupted power supply (UPS).

Remove the connector plug, push the release button, and then unplug the connector plug.

2.3. Connect the pH and conductivity electrodes

The following graphics demonstrate sockets available for pH and conductivity measurement. Some models have fewer sockets than others, depending on the measurement modes the meter/analyzer support.

a-AB33M1

There are two sockets for pH electrode, the “pH/mV” socket and the “Temp” Socket. One socket for conductivity electrode, the “COND” socket.

Temp



Figure 2-2 a-AB33M1 sockets

a-AB33PH

There are two sockets for pH electrode, the “pH/mV” socket and the “Temp” Socket.

Note: For 2-in-1 pH electrode, you only need to connect it to the pH/mV socket, while for other electrodes, you need to connect to both the pH/mV and Temp socket.

a-AB33EC

There is one socket “COND” for conductivity electrode.


Figure 2-4 a-AB33EC sockets

2.4. Touch screen

AquaSearcher™ series desktop meter/analyzer touch screen, including menu-related function keys, enables it to operate the meter/analyzer intuitively and efficiently, and easily navigate to the menu of measurement, setting, data log and calibration log.


Figure 2-5 Touch screen

1. **Date & Time**
2. **Clock icon:** when the reading type is set to timing function, it is always on.
3. **Computer icon:** it indicates that the meter/analyzer has been connected the computer, and if it disappears, it means that the connection is disconnected.
4. **Auto-Stop icon ((^)):** when the reading type is set to Auto-Stop, the icon will blink when reading is stable; when reading is stable and the display is locked, the icon is stable.
Continuous reading icon((c)): when the reading type is set to self-continuous reading, the icon will blink when the reading is stabilizing.
5. **Measurement icon:** blinks when measurement is in progress, and stops blinking when measurement is finished.
6. **Printer icon:** touch to print.
7. **Calibration icon:** always on indicates that a calibration has been successfully completed; blinking indicates that the calibration alarm is set and not calibrated; touch to calibrate.
8. **Menu icon:** touch to enter the menu settings.
9. **pH electrode status (a-AB33M1, a-AB33PH)**

Slope: 95-105%
and offset: (0-15) mV.

Slope 90-95%
or offset: (15-35) mV.

Slope: less than 90%
or offset: (35-60) mV.

Electrode condition is good.

Electrode condition is acceptable.

Electrode condition is not good or needs cleaning.

Conductivity electrode condition (a-AB33M1, a-AB33EC)

Calibration result: within $\pm 20\%$ of the set cell K value.

Electrode condition is good.

Calibration result: out of $\pm 20\%$ of the set cell K value.

Electrode needs cleaning.

10. First measurement channel (a-AB33M1)

Note: a-AB33M1 supports dual channel measurement, other models only support one measurement channel.

11. **Temperature:** touch to switch between automatic temperature compensation (ATC)/ manual temperature compensation (MTC).

12. **Measurement mode area:** touch to switch between different pH and conductivity measurement modes.

pH, mV, ORP (a-AB33M1, a-AB33PH);

EC, TDS, SAL (Salinity), RDS (Resistivity) (a-AB33M1, a-AB33EC);

13. Second measurement channel (a-AB33M1)

Note: a-AB33M1 supports dual channel measurement, other models only support one measurement channel.

14. **Beaker icon:** touch to recall measurement data.

15. **User ID icon:** touch and press the or button to select user by user ID. The main interface is displayed at the lower part of the screen in measurement state: Summary of the latest calibration.

Touch the **Menu** icon  to enter the following menu: touch the F1, F2 and F3 function keys to execute the functions displayed above the keys displayed on the display screen.

16. **F1 Setup:** touch to enter the parameter setup menu.

17. **F2 i-Steward:** touch to enter electrode condition menu.

18. **F3 Settings:** touch to enter meter's setup menu.

Note: when you enter the three menus introduced above, touch the related sub-menu items to continue settings.

2.5. Buttons

Up button:

- When setting parameters, press to increase the setting value or to display the previous option.

- b) When in setup menu, press to scroll up to the previous page of menu selections.
- c) When in calibration and PH electrode check system (PH ECS), press to repeat the previous step.

**Home button:**

- a) Press for about three seconds to turn on the meter.
- b) Press to start a new measurement.
- c) Press to exit from the current mode or menu and return to the measurement mode.
- d) Press for about three seconds to turn off the meter.

**Down button**

- a) When setting parameters, press to decrease the setting value or to display the next option.
- b) When in setup menu, press to scroll down to the next page of menu selections.
- c) When in calibration and PH electrode check system (PH ECS), press to move to the next step.
- d) Change the measurement channel displayed.

2.6. Turn on and off the meter/analyzer

Turn on the meter/analyzer

1. Touch  Home button for 3 seconds and the screen will light up.
2. The meter/analyzer is initialized and the system is tested, and it will take some time for this process.
3. Touch  or  button to select user ID, and touch  user ID.
 - a) If no password is set for the user ID, enter the main interface directly.
 - b) If a password has been set for the user ID, enter the 4-digit password to enter the main interface.

Turn off the meter/analyzer

1. Touch  Home button for 3 seconds and the screen will go out. When the screen goes out, the meter/analyzer is in standby state.

Note: In standby mode, the control circuit of the  Home button is energized. Other parts of the meter/analyzer will not be powered on.

3. User management

User management helps to manage measurement results and parameter settings by different users.

The meter/analyzer provides 10 users. User ID 0 refers to the admin account, which has the function of clearing records, but it is mandatory to export records before clearing log.

User management cannot be accessed during measurement.

3.1. User management

To switch between different users:

1. Touch  on the screen.
2. Press the  or  to select between different user IDs. Touch  when you see the ID number you want to confirm your setting.

4. Records management

 The meter/analyzer uses the menu to access the measurement records. All the information in the measurement record is saved in the memory of the meter/analyzer, so the saved information will not be lost even if it is not connected to the power supply. The meter/analyzer can record 1000 sets of measurement results. Each record contains the measurement results displayed actively, depending on the meter/analyzer model and display configuration, with a date and time stamp.

Navigation: screen icon .

- 1. Touch  or  button to view each record. If total number of data records exceeds 10, 100, and 500, the screen displays the quick navigation, and it can scroll to the 10th, 100th and 500th records.
- 2. If you connect the meter/analyzer to a printer or computer, touch  on the screen to print the record.
- 3. Touch  again to exit and return to the measurement screen.

Total number of measurement and calibration records supported by AB33 series meter/analyzer.

Table 4-1 Records management

Model	Maximum Measurement Records	Maximum Calibration Records
AB33PH	1000	1 pH calibration result and 1 ORP calibration result
AB33EC	1000	1 conductivity calibration result
AB33M1	1000	1 pH calibration result, 1 ORP calibration result and 1 conductivity calibration result

Note:

- The calibration records are separate and not included in the 1000 sets of measurement records.
- Continuous measurement mode, the measured data results will not be recorded.
- If the record is full, the oldest record will be automatically overwritten.

5. Setup menu

Touch the **Menu** icon  to enter the meter/analyser’s setup menu. For the menu structure, refer to the following:

Menu					
	Setup			i-Steward	Settings
Mode	pH	Conductivity	Temperature	Cal. due alarm	Language
Read type	Resolution	Cell K	Temperature Unit	Export Log	RS232
Stability	Buffer Group	Std. Solution	Enter temp. value	Log Clear	Date/time
Average		Ref. Temp.	Temp. calibration	Self-Test	Brightness
		Temp. Coeff.		Reset	Keypad Beep
		TDS Factor			Auto Shut Off
					About Meter

5.1. Setup

Enter this menu to set the menu of related measurement parameters, such as setup mode, pH channel measurement parameters, conductivity channel measurement parameters and temperature setting. It is set as bold by default

Table 5-1 Setup menu

Menu	Sub-Menu	Sub-Menu Options	Option
Setup	Mode	Read type	Auto-Stop Continuous reading
		Stability	SmartStability Fast Medium
		Smart Avg.	OFF On
	pH	Resolution	0.1 0.01

Menu	Sub-Menu	Sub-Menu Options	Option
	(a-AB33M1, a-AB33PH)	Buffer Group	USA JJG DIN EU Merck
	Conductivity (a-AB33M1, a-AB33EC)	Cell K	0.01~10.0 (the default is 0.1 (a-AB33EC) and 1.75 (a-AB33M1))
		Std. Solution	10 µS/cm 84 µS/cm (a-AB33EC) 146.5 µS/cm 500 µS/cm 1413 µS/cm (a-AB33M1) 12.88 mS/cm
		Ref. Temp.	20°C (68°F) 25°C (77°F)
		Temp. Coeff.	0.0~10.0%/°C (the default is 2.0%/°C)
		TDS Factor	0.01~5.0 (the default is 0.50)
	Temperature	Temperature Unit	°C °F
		Enter temp. value	5.0-60.0°C (41.0-140°F) (the default is 25.0°C/77.0°F)

5.1.1 Mode

5.1.1.1 Read type

Navigation: Menu  → Setup → Mode → Read type

Select the Read type to determine the setting of the measurement end mode. The Read type determines how the meter/analyzer displays, records logs and exports measurements.

Auto-Stop = the meter/analyzer determines when the reading is stable and reaches the endpoint automatically.

Continuous = the meter/analyzer keeps measuring and users manually end the measurement.

Auto End Point Mode

1. Touch the  button to start measurement.  icon start blinking during measurement.
2.   icon blinking twice indicates that the reading has reached the endpoint. The results will be automatically stored. **Note:**

- Touch the  icon to print the result if you connect the meter to a printer or a computer.
- Touch  button to start the new measurement.

Continuous Reading Mode

1. Touch the  button to start measurement.  icons and  icon start blinking during measurement.

The blinking icon  indicates that if you connect the PC to the meter/analyzer, the measurement results will be continuously output to the PC, and the results will not be stored in the meter/analyzer.

Note: During the output, the frequency of the pH channel is 2Hz, while the frequency of the conductivity channel is 1Hz.

2. Touch  button to end the measurement.

5.1.1.2 Stability

Navigation: Touch the menu icon  → Setup → Mode → Stability

The algorithm of stability to reading fluctuation determines that the meter/analyzer measurement is recognized as stable. After selecting “Smart Stability”, the measurement will be ended in combination with the algorithm and measurement resolution, and the results will be displayed.

SmartStabilit
y = the meter/analyzer determines stabilization time automatically.

Fast = faster stabilization time with less stability.

Medium = normal stabilization time with normal stability

The table shows the relationship between estimated millivolts and time under ideal measurement conditions. This algorithm is only a theoretical value, and the actual measurement results will be different depending on the measurement conditions.

Table 5-2 Stability algorithm

Select	Algorithm
Smart Stability (with 0.1 resolution or 1 significant digit)	When the difference between the maximum value and the minimum value of the measurement result within 4 seconds is less than 0.5pH or 0.6mV, the meter/analyzer will consider the result as stable.
Smart Stability (with 0.01 resolution or 2 significant digit)	When the difference between the maximum value and the minimum value of the measurement result within 6 seconds is less than 0.08pH or 0.15mV, the meter/analyzer will consider the result as stable.
Fast stability	When the difference between the maximum value and the minimum value of the measurement result

	within 4 seconds is less than 0.5pH or 0.6mV, the meter/analyzer will consider the result as stable.
Medium stability	When the difference between the maximum value and the minimum value of the measurement result within 6 seconds is less than 0.08pH or 0.15mV, the meter/analyzer will consider the result as stable.

Note:please enter the following PH resolution menu to set resolution mentioned in the first three selections in the table.

5.1.1.3 **Smart Avg.**

Navigation: Touch the menu icon  → Setup → Mode → Smart Avg.

The average measurement result will be calculated by the meter/analyzer automatically. This selection helps to achieve stable measurements more rapidly.

- On = enabled
- OFF = disabled The measurements will not be averaged after meeting the stabilization criteria.

5.1.2 **pH (a-AB33M1, a-AB33PH)**

a-AB33M1 and a-AB33PH meter/analyzer customize the resolution and buffer group for the pH channel.

5.1.2.1 **pH resolution**

Navigation: Touch the menu icon  → Setup → pH → resolution setting, resolution of pH measurement.

- 0.1 = 1 decimal place
- 0.01 = 2 decimal place

5.1.2.2 **Buffer Group**

This section introduces the setting of pH preset buffer group in a-AB33M1 and a-AB33PH meter/analyzers.

During the pH calibration, the meter/analyzer uses the selected pH buffer group and the original mV value of the pH electrode in the buffer to identify and display the buffer value of the measured temperature. The original mV value of the pH electrode in the buffer solution must be within one pH unit (about 59 mV) of the theoretical MV value of the buffer solution of the meter/analyzer to identify the buffer solution automatically.

Before calibration, the preset buffer group must be defined, and the preset buffer group and the standard solution in the group cannot be added, modified or deleted.

Navigation: Touch the menu icon  → Setup → pH → Buffer group

The following parameters of buffer and standard solution have been stored in the meter/analyzer, and the buffer is automatically identified. Please refer to the appendix.

Table 5-3 Buffer group

pH Buffer Group	Reference temperature (Ref. Temp.)	pH Value				
		1.68	4.01	7	10.01	/
USA standard	25°C	1.68	4.01	6.86	9.18	12.46
JJG standard	25°C	2	4.01	7	9.21	11
EU standard	20°C	2	4	7	9	12
MERCK standard	25°C	1.68	4.008	6.865	9.184	12.454
DIN (19266) standard	25°C	1.68	4.008	6.865	9.184	12.454

5.1.3 Conductivity channel setting (a-AB33M1, a-AB33EC)

5.1.3.1 Cell K

The Cell K value is usually recorded in the factory inspection report or marked on the electrode's cable, which is about 0.01~10.0 cm⁻¹. The default value of a-AB33EC model is 0.1 cm⁻¹, and that of a-AB33M1 model is 1.75 cm⁻¹.

Navigation: Touch the menu icon  → Setup → Conductivity → Cell K compatible conductivity Cell K value:

Table 5-4 Cell K

Electrode model	Cell K value
STCON8	0.08-0.1
STCON7	0.08-0.1
STCON5	0.8-1.2
STCON3	1.5-2.0

After calibration, if the result is within ±20% of the Cell K value you set, the screen will display , indicating that the electrode is in good condition. If the result is out of range, the screen will display , indicating that the electrode needs to be cleaned.

5.1.3.2 Standard solution (Std.Solution)

Navigation: Touch the menu icon  → Setup → Conductivity→ Std. Solution

The six predefined standard solutions are 10 μS/cm, 84uS/cm, 146.5 μS/cm, 500 μS/cm, 1413uS/cm and 12.88 mS/cm respectively.

5.1.3.3 Reference temperature (Ref. Temp.)

Navigation: Touch the menu icon  → Setup → Conductivity→ Ref. Temp.

Select the range: 20°C (68°F)

25°C (77°F)

5.1.3.4 Temperature coefficient (Temp. Coeff.)

Navigation: Touch the menu icon  → Setup → Conductivity → Temp. Coeff.

Input range: 0.0-10.0%/°C (2.0%/°C by default)

Note: If you set the temperature correction coefficient value to 0.0 %/°C, there is no temperature compensation for conductivity measurement. The meter/analyzer will display the true conductivity value at the current uncompensated temperature.

5.1.3.5 TDS factor

Navigation: Touch the menu icon  → Setup → Conductivity → TDS Factor

Input range: 0.01-5.0 (0.50 by default)

Touch  or  button to adjust the temperature value. Touch the values on the screen to confirm the compensation temperature settings.

5.1.4 Temperature

The a-AB33 series meter/analyzer provides a temperature menu, which allows you to select °C or °F as the temperature unit, customize the temperature compensation setting for each channel by manually entering the temperature value, and perform temperature correction for ATC probe, conductivity sensor or dissolved oxygen sensor with built-in temperature.

The meter/analyzer is connected to the temperature electrode, and the automatic temperature compensation (ATC) and sample temperature will be displayed on the screen. You can choose to skip the manual temperature compensation (MTC) setting.

If the temperature electrode is not detected or not used by the meter/analyzer, the meter/analyzer will automatically switch to the manual temperature compensation (MTC) mode and display MTC on the screen.

5.1.4.1 Temperature unit

Navigation: Touch the menu icon  → Setup → Mode → Temperature → Temp. unit

°C = °C

°F = °F

Calculation formula of conversion from °F to °C: $^{\circ}\text{C} = 5/9 (^{\circ}\text{F} - 32)$.

5.1.4.2 Enter temperature value (Enter temp. value)

A-AB33 series meter/analyzer is configured with manual temperature entering. It is set according to the actual sample temperature, and a temperature compensation is provided. Navigation: Touch the menu icon  → Setup → Mode → Temperature → Enter temp. value

Input range: 5.0~60.0°C (41.0~140°F) (25.0°C/77.0°F by default)

Touch  or  button to adjust the temperature value. Touch the values on the screen to confirm the compensation temperature settings.

5.2. i-Steward™ menu

Enter this menu to set electrode condition parameters.

Table 5-5 i-Steward menu

Menu	Sub-Menu	Option
i-Steward	Cal. due alarm	OFF 1 day 7 day 14 day
	pH ESC	Enter the pH ECS process, please complete the operation according to the instructions.
	Export Log	Result Log Cal. Log
	Log Clear (ID0 only)	Result Log Cal. Log
	Self-Test	/
	Reset	Factory reset Parameter Reset

5.2.1 Calibration due alarm (Cal. due alarm)

A-AB33 series meter/analyzer provides Cal. due alarm to remind customers to check electrode performance. After turning on, for the reminding of calibration expiry, the calibration interval can be set by day. If no calibration is performed within the calibration period set, an alarm will be triggered, and the calibration icon  will blink, until it is always on after completing the calibration. Enter this menu to enable or disable system the Cal. due alarm. Navigation: Menu  → i-Steward™ → Cal. due alarm

- OFF = disabled
- 1 day = the meter/analyzer will warn if no calibration within 1 day after the last calibration.
- 7 day = the meter/analyzer will warn if no calibration within 7 day after the last calibration.
- 14 day = the meter/analyzer will warn if no calibration within 14 day after the last calibration.

5.2.2 pH Electrode Check System (pH ECS)


Enter this menu to start the pH ECS process, please complete the operation according to the instructions.

Navigation: Menu  → i-Steward™ → pH ESC


The process is as follows:

1. Take out the optional electrode inspection tester connector and remove the BNC short circuit cap on the connector.


Note: You can contact an authorized Ohaus dealer to purchase the connector of the electrode inspection tester.

2. Connect the electrode inspection tester to the pH/mV socket of the meter/analyzer as shown in the figure below.

3. Prepare the electrodes as displayed on the screen, and then press the button  to continue when finished.
4. Wait for the reading to stabilize. When you see the screen displaying “parallel 100MΩ resistors in buffer solution with pH4.01”, close the BNC short-circuit cap on the connector as shown in the figure below. Press the button  to continue when finished.

5. The display screen displays the results. The meter/analyzer will display whether your electrodes are in good condition.

5.2.3 Export log

a-AB33 series meter/analyzer provides logs for export to PC. Enter this menu to export all Result Logs and Cal. Log.

Please use SPDC serial port and Ethernet data acquisition tool (PC-end software) to connect the meter/analyzer for receiving and export, and generate a PDF report.

Navigation: Menu  → i-Steward™ → Export Log

Select the specific categories of records that need to be exported, namely Result Log and Cal. Log.

Note: Select the meter/analyzer type as AquaSearcher.

For better compliance, export the document format as PDF.

5.2.4 Log clear (ID0 only)

A-AB33 series meter/analyzer can clear all records, including Result Logs and Cal. Logs. This function is only visible to the administrators (user ID0).

Navigation: Menu  → i-Steward™ → Log Clear

1. Select the specific categories of records that need to be cleared, namely Result Log and Cal. Log.
2. Touch  button to confirm deletion, or touch  button to cancel deletion. Wait for the display process to be completed. **Note:** When the screen displays “Cannot clear before export.”, operate with reference to **5.2.3 Export log**
3. After the clearing is completed, touch  button or  button to return to the Setup menu.

5.2.5 Self-test

A-AB33 series meter/analyzer provides hardware self-inspection function. Enter this menu to test if the touch screen, and the buttons function normally.

Navigation: Menu  → i-Steward™ → Self-Test

1. If you use a-AB33M1 and a-AB33PH, please disconnect all electrodes and connect the supplied BNC short-circuit cap to the pH/mV socket of the meter/analyzer (if applicable). Touch the button  to continue when finished. If you are use other models, please skip this step.
2. Wait for the reading to stabilize and then enter the testing the screen step. Check the display to see if all the icons are displayed normally.
3. Press the  button to enter the next step to test the touch screen. Touch all the icons on the screen one by one. Their disappearances from the screen indicate that the touch screen functions normally.
4. Touch the  or  button one by one according to reminder on the screen. If all the above steps have been passed, the meter/analyzer will show a reminder that the self-test has passed.
5. Press the  or  button to return to the main menu.

If your meter/analyzer does not pass the self-test and show errors. You can contact an OHAUS technical support agent for help.

5.2.6 Reset

A-AB33 series meter/analyzers provide reset function, which can restore factory settings and parameter settings.

5.2.6.1 Factory reset

Warning: Performing factory reset will delete all measurement data logs and calibration data logs of the meter/analyzer and reset all meter/analyzer setting parameters to factory default settings.

Enter this sub-menu to reset the meter/analyzer to its factory default settings. This function is only visible to the administrators (user ID0). Touch the following buttons in turn according to the navigation:

Navigation: menu  → i-Steward™ → Reset → Factory reset

Note: Before factory reset, when the screen displays “Cannot clear before export.”, the meter/analyzer is forced to execute the export record first, operate with reference to **5.2.3 Export log**.

- 1. Press the  button to confirm, or press the  button to cancel.

5.2.6.2 Parameter reset

Warning: Performing parameter reset will reset to factory default settings; however, the calibration log and data log will remain. Enter this sub-menu to reset all the parameters under the Setup menu. Touch the following buttons in turn according to the navigation:

Navigation: Menu  → i-Steward™ → Reset → Parameter reset

Press the  button to confirm, or press the  button to cancel.

5.3. Settings menu

Enter this menu to set the meter/analyzer’s settings. The default settings are in bold.

Table 5-6 Settings menu

Menu	Sub-Menu	Option	
Settings	Language	English	
		Espanol	
		Francais	
		中文	
		Türkçe	
		한국인	
		Русск	
		Portuguese	
	RS232	Export to	PC , printer
		Print Length	Complete Length, Simple Length
		Baud rate	38400, 19200, 9600, 4800, 2400, 1200
		Parity	7 Even
			7 Odd
			7 None
		Stop Bits	8 None
			1 bit
		Handshake signal	2 bit
			None
		Xon/Xoff handshake	

	Date/time	Time format	24H Clock 24H Clock
		Date format	MM/DD/YY DD/MM/YY YY/MM/DD
		Adjust Date & Time	/
	Brightness	High Medium Low	
	Keypad Beep	OFF On	
	Auto Shut Off	OFF 15 min 30 min	
	About Meter	/	

5.3.1 Language

Select the language to be used in the meter/analyzer interface, and set the language of menu and information display.

Navigation: Menu  → Settings → Language

- English
- Espanol
- Francais
- 中文
- Türkçe
- 한국인
- pyccк
- Portuguese

5.3.2 RS232

To transmit the measurement results to a PC or printer, update the export settings, including output format, baud rate, parity check bit, etc. according to the settings of external devices in the sub-menu.

Navigation: Menu  → Settings → RS232

Export to

Set the position of output data:

- Printer = export to printer, with transmission format of the printer format.
- PC = export to a computer, with transmission format of CSV.

Print Length

- Complete Length = Print header and test record data.
- Simple Length = Print test record data

Baud rate

Set Baud rate (bit per second)

1200	= 1200 bps
2400	= 2400 bps
4800	= 4800 bps
9600	= 9600 bps
19200	= 19200 bps
38400	= 38400 bps

Parity

Set data bit and parity check bit

7 Even	= 7 data bits, even parity
7 Odd	= 7 data bits, odd parity
7 None	= 7 data bits, no parity
8 None	= 8 data bits, no parity

Stop Bits

Set the number of stop bits.

1 bit	= 1 stop bits
2 bit	= 2 stop bits

Handshake signal

Set the flow control method.

None	= no handshaking
Xon/Xoff	= XON/XOFF software handshaking

5.3.3 Date/time

A-AB33 series meter/analyzer can provide relevant settings for date and time. This sub-menu enters a display format for the date and time and set the current date and time.

Navigation: Menu  → Settings → Date/time

Time format

Set the time format

24H Clock	= 12-hour format, displayed as AM or PM
24H Clock	= 24-hour format

Date format

Set date format

MM/DD/YY	= MM/DD/YY
DD/MM/YY	= DD/MM/YY
YY/MM/DD	= YY/MM/DD

Adjust Date & Time

Enter this menu to adjust date and time.

1. The first two digits of the date at the upper left corner of the screen start to blink.
2. Touch the  or the  button to adjust to the current date and time. When finished, touch Next on the screen, and all digits have been adjusted.
3. Touch the menu , and the meter/analyzer returns to the main menu.

5.3.4 Brightness

Enter this menu to increase or decrease the display brightness under different lighting conditions.

Navigation: Menu  → Settings → Brightness

High	= high brightness
Medium	= medium brightness
Low	= low brightness

5.3.5 Keypad beep

Enter this menu to enable or disable the beep sound when pressing buttons.

Navigation: Menu  → Settings → Keypad Beep

On	= enabled beep sound
OFF	= disabled beep sound

5.3.6 Auto shut off

Enter this menu to set up automatic shutdown time. The meter/analyzer will shut down if no operation during the time period you set.

Navigation: Menu  → Settings → Auto Shut Off

OFF	
15	= the meter/analyzer will shut down after 15 minutes of no operation.
30	= the meter/analyzer will shut down after 30 minutes of no operation.

Note:during auto or continuous measuring modes, the meter/analyzer will not perform this operation.

5.3.7 About meter

Enter this menu to check the meter's model, serial number and software version.

Navigation: Menu  → Settings → About Meter

6. pH operation (a-AB33M1, a-AB33PH)

This chapter is for a-AB33M1 and a-AB33PH models which support pH measurement. Before you start direct calibration, you must set the following parameters:

- Refer to Chapter 5 Mode in meter/analyzer setting menu
- Resolution and buffer group in pH parameter setting menu
- Temperature setting and data transmission in the Settings menu
- Data log
- Date & Time
- Language
- Sound
- Brightness
- Auto Shut Off, etc.

6.1. pH calibration process

a-AB33M1 and a-AB33PH meter/analyzer models support a maximum 5-point pH buffers calibration, and at least two-point calibration is recommended according to the pH or range of the sample.

Standard procedure of pH measurement is as follows:

1. pH electrode preparation: Prepare pH electrodes as instructed in the electrode manual, and connect all the electrodes to be used to the corresponding meter/analyzer inputs. pH electrode should be rinsed with pure water before and after use. Check if the electrode is physically damaged.
2. In buffer preparation, be sure to use fresh pH buffer and choose a buffer containing the sample pH and one to three pH units away;
3. In pH measurement, confirm that the meter/analyzer measurement mode has been switched to pH channel;
4. Record the measurement results or print;
5. Rinse the pH electrode and properly store;
6. The pH electrode should be stored in the storage bottle; the solution in the bottle is 3M KCl solution.

Warning Do not operate the equipment in hazardous environments. The equipment is not explosion protected.

Warning When using chemicals and solvents, comply with the instructions of the chemical producer and the general lab safety rules.

6.1.1 pH electrode direct calibration process

Select a pH resolution and a predefined buffer group before performing the calibration. For details, please refer to **5.1.2 pH** above. Before starting calibration, prepare the electrode as instructed in the electrode manual. Connect the electrode to the meter/analyzer.

1. Confirm that the measurement channel has been switched to pH channel, and touch  **Cal** icon on the screen to start calibration.
2. Rinse the pH electrode and all other electrodes with pure water as instructed on the screen. Wipe off the excess water with a lint-free cloth and put it in the calibration buffer.
3. When the electrodes and buffer are ready, touch the  button to start the calibration. The  icon on the screen starts to blink during the process. **Note:** The calibration starts from Auto End Point Mode.
4. Wait until the pH value on the meter/analyzer is stable, and the meter/analyzer reaches the endpoint automatically. The screen will display the calibration point pH value (for example, 7.00) with the temperature. Then execute one of the following operations:
 - Touch **Accept result** to accept the pH buffer value identified by the current automatic buffer.
 - Touch **Re-Calibration** to re-calibrate this point, and repeat the previous step.
 - Touch  button to cancel the calibration process, without saving the results, and exit the process.
5. Touch **Accept results**, and then execute one of the following operations:
 - Touch "Finish" to complete the calibration, and the calibration result will be displayed.
 - Touch  to continue to the next pH buffer and repeat steps 2 to 5. If five buffers are used, the calibration will be saved and ended when the fifth buffer is accepted.
The calibration summary (including slope, offset and calibrated buffer point) will be displayed at the bottom of the screen, and the data will be automatically saved in the calibration log. Touch  button to exit the calibration and return to the measurement interface, or if you connect the meter to a printer or PC, touch  icon to print the results.

6.1.2 ORP electrode (relative mV mode) calibration process

In **ORP** mode, use a standard solution for single point calibration, and perform ORP calibration according to the on-screen instructions. Before starting calibration, prepare the electrode as instructed in the electrode manual. ORP standard solution is non-hazardous, stable and will not change with time, and fresh ORP standard solution will be used for validation.

Connect the ORP electrode to the meter/analyzer for ORP calibration:

1. Touch to select ORP measurement mode, confirm that the unit has been set to relative mV unit, and touch  **Cal** icon on the screen.
2. Rinse with pure water as instructed on the screen. Wipe off the excess water with a lint-free cloth and put it in the standard solution.

3. When the electrodes and buffer are ready, touch the  button to start the calibration. The  icon on the screen starts to blink during the process. **Note:** The calibration starts from Auto End Point Mode.
4. Wait until the relative mV value on the meter/analyzer is stable, and the meter/analyzer reaches the endpoint automatically. Touch "Finish" to complete the calibration, and the calibration result will be displayed. (For example, 223mV, ± 3 mV, 230mV and 25°C).
 - Touch **Recalibration** if you want to re-calibrate and repeat the previous steps.
 - Touch  button to cancel the calibration process and return to the measurement interface, without saving the results, and exit the process.
 - Touch **Finish** to complete the calibration, and the calibration result will be displayed. The calibration summary (including slope, offset and calibrated buffer point) will be displayed at the bottom of the screen, and the data will be automatically saved in the calibration log.
 - If you connect the meter/analyzer to a printer or PC, touch  icon to print the results.

6.2. pH measurement process

Prepare the electrode to be used as instructed in the electrode manual. Connect all the electrodes to be used to the corresponding meter/analyzer inputs. Make sure that the electrode has been calibrated recently and is operating normally.

Make sure that the meter/analyzer measurement mode is set to the required measurement parameters.

6.2.1 pH measurement

1. Wash the electrode with distilled water or a suitable solution, wipe off the excess water gently with a lint-free cloth, then put it in the sample to be tested, and then wait for 30 seconds.
2. Touch  button to start the pH measurement, and during measurement,  icon will blink.
3. After the measurement is stable or reaches the set time ( icon is displayed on the display screen, and  icon stops blinking), the meter/analyzer will display the pH value and temperature, and record all relevant data.
 - Auto-Stop mode: when the measurement is stable, it will be locked to the display screen. Touch the  icon to print the result if you connect the meter to a printer or a computer.

- Continuous endpoint mode: after the measurement is stable, touch  button to end the measurement manually. Connect the meter/analyzer to a printer or PC and transmit the measurement data in real time.
 - Timing measurement mode: record the measurement at predetermined time intervals. Touch the  icon to print the result if you connect the meter to a printer or a computer.
4. Take out the electrode from the sample, wash the electrode with distilled water or a suitable solution, gently wipe off the excess water with a lint-free cloth, and then put it in the next sample to be tested.
 5. Repeat steps 2 to 4 for all samples. After measuring all samples, keep the electrodes as instructed in the electrode manual.

6.2.2 mV and ORP measurement

- 1 Touch the mode area on the screen to switch between mV and ORP measurement.
- 2 Follow the same measurement procedure as pH measurement.

6.2.3 Temperature measurement

For a better measurement accuracy, we suggest using a three-in-one or separate temperature electrode. Use linear temperature compensation.

- If a temperature electrode is used, ATC and sample temperature will be displayed.
- If the meter/analyzer does not detect the temperature electrode, it will automatically switch to manual temperature compensation mode and MTC will appear. MTC temperature should be set before temperature compensation.

Note: a-AB33PH and a-AB33M1 models support the use of MTC mode, even if temperature electrodes are used. If MTC mode is selected, MTC will continue to be displayed.

Note: The meter/analyzer accepts NTC 30k Ω temperature electrode. ATC (automatic temperature compensation) or MTC (manual temperature compensation): Temperature compensation only corrects the change of electrode output, not the change of actual solution. This means that the meter/analyzer will correct the signal (mV) from the electrode according to the actual temperature to obtain a more accurate pH value.

7. Conductivity operation (a-AB33M1 and a-AB33EC)

This section applies to the a-AB33M1 and a-AB33EC models that support conductivity measurement.

7.1. Conductivity calibration process

- 1 Before calibration, set calibration standard solution and Cell K. For details, please refer to **5.1.3.1 Cell K** and **5.1.3.2 Standard solution** above. Touch the mode area one or more times to select **EC**.
- 2 Touch  icon on the screen.
- 3 If you use a-AB33M1 model, please touch **Conductivity**. If you are use other models, please skip this step.
- 4 Rinse the electrode as instructed on the screen, put it in the standard solution, stir for 5 seconds and wait for 30 seconds.
- 5 Press the  button when you finish the previous step to start calibration. The  icon on the screen starts to blink during the process.

Note:

- Calibration starts with the Auto End Point Mode.
 - Pressing the  button during the calibration can cancel the process.
- 6 The meter/analyzer reaches the endpoint automatically. The screen will display the calibration point EC value with the temperature.
 - Touch **Finish** to complete the calibration, and the calibration result will be shown.
 - Touch **Recalibration** if you want to re-calibrate and repeat the previous steps.

Note:

- Pressing the  button can cancel the calibration and return to measurement screen.
- Touch the  icon to print the result if you connect the meter to a printer or a computer.

7.2. Conductivity measurement

Prepare the electrode to be used as instructed in the electrode manual. Connect all the electrodes to be used to the corresponding meter/analyzer inputs. Make sure that the electrode has been calibrated recently and is operating normally.

Make sure that the meter/analyzer measurement mode is set to the required measurement parameters.

- 1 Place the electrode in the sample, stir for 5 seconds and then wait for 30 seconds
- 2 Press the  button to start the conductivity measurement, the  icon starts to blink during the measurement.
- 3 When the meter reaches the end point (the  icon appears on the display and the  icon stops blinking), it will display the conductivity value and the temperature.

Note: in Continuous End Point Mode, you need to press the  button to manually end the measurement.

7.3. TDS, SAL (Salinity) and RES (resistivity) measurement

- 1 Touch the mode area on the screen one or several times to switch between **TDS**, **SAL** and **RES** modes.
- 2 Follow the same procedure as for conductivity measurement.

8. Print

8.1. Format

Note: in the following tables, number and item columns are only for illustration to explain the printing content and will not be printed out in the final result.

Table 8-1 pH calibration

No.	Item		Printing Content
1	Meter Model		AB33PH/AB33M1
2	Serial Number		XXXXXXXXXX
3	Software Revision		X.XX
4	Channel		pH
5	Mode		pH
6	Date & Time		September 25, 2020, 14:34
7	Calibration points		1-5
8	Repeat ed for each point-to- point	Cal Point Index	1
9		pH Value	X.XX
10		pH Unit	pH
11		mV Value	XX.X
12		mV Unit	mV
13		Temperature Value	XX.X
14		Temperature Unit	C°, F°
15	Slope Value		XXX.X
16	Slope Unit		%
17	Offset Value		XXX.X
18	Offset Unit		mV

Table 8-2 ORP calibration

No.	Item		Printing Content
1	Meter Model		AB33PH/AB33M1
2	Serial Number		XXXXXXXXXX
3	Software Revision		X.XX
4	Channel		pH
5	Mode		ORP
6	Date & Time		September 25, 2020, 14:34
7	ORP mV Value		XXX.X
8	ORP mV Unit		mv
9	Offset mV Value		XX.X
10	Offset mV Unit		mv
11	Reference mV Value		XXX.X
12	Reference mV Unit		Rmv
13	Temperature Value		25
14	Temperature Unit		C°, F°
15	Temperature Type		ATC/MTC
16	Calibration ID		#1

Table 8-3 Conductivity calibration

No.	Item	Printing Content
1	Meter Model	AB33EC/AB33M1
2	Serial Number	XXXXXXXXXX
3	Software Revision	X.XX
4	Channel	Conductivity
5	Mode	Conductivity
6	Date & Time	September 25, 2020, 14:34
7	Conductivity Value	XX.X
8	Conductivity Unit	uS/cm
9	Conductivity Value	XX.X
10	Conductivity Unit	uS
11	Temperature Value	XX.X
12	Temperature Unit	C°, F°
13	Temperature Type	ATC/MTC
14	Calibration Type	Manual
15	Cell K	X.X
16	Calibration ID	1

Table 8-4 pH measurement

No.	Mode		Printing Content
	pH	mV	
1	Meter Model		AB33PH/AB33M1
2	Serial Number		XXXXXXXXXX
3	Software Revision		X.XX
4	User ID		0-9
5	Date & Time		September 25, 2020, 14:34
6	Sample ID		000-999
7	Channel		pH
8	Mode		pH, mV
9	pH Value	/	XXX.X
10	pH Unit	/	pH
11	mV Value		XXX.X
12	mV Unit		mV
13	Temperature Value		XX.X
14	Temperature Unit		°C, °F
15	Temperature Type		ATC, MTC
16	Slope Value		XXX.X
17	Slope Unit		%
18	Offset Value		XXX.X
19	Offset Unit		mV

Table 8-5 ORP measurement

No.	Item	Printing Content
1	Meter Model	AB33PH/AB33M1
2	Serial Number	XXXXXXXXXX
3	Software Revision	X.XX

4	User ID	0-9
5	Date & Time	September 25, 2020, 14:34
6	Sample ID	000 -999
7	Channel	pH
8	Mode	ORP
9	ORP Value	XXX.X
10	ORP Unit	RmV
11	Offset mV Value	XXX.X
12	Offset mV Unit	mV
13	Temperature Value	XX.X
14	Temperature Unit	C°, F°
15	Temperature Type	ATC, MTC

Table 8-6 Conductivity measurement

No.	Mode				Printing Content
	Conductivity	TDS	SALT	RES	
1	Meter Model				AB33EC/AB33M1
2	Serial Number				XXXXXXXXXX
3	Software Revision				X.XX
4	User ID				0-9
5	Date & Time				September 25, 2020, 14:34
6	Sample ID				000-999
7	Channel				Conductivity
8	Mode				Conductivity, TDS, SALT, RES
9	Conductivity Value	TDS Value	SALT Value	RES Value	XXX.X, XX.XX, X.XXX
10	Conductivity Unit	TDS Unit	SALT Unit	RES Unit	uS/cm, mg/L, psu, Ω.cm
11	Conductivity Value				XXX.X, XX.XX, X.XXX
12	Conductivity Unit				S
13	Temperature Value				XX.X
14	Temperature Unit				C°, F°
15	Temperature Type				ATC, MTC
16	Temp. Coeff.	Temp. Coeff.	/	Temp. Coeff.	XX.X
17	Temp. Coeff. Unit	Temp. Coeff. Unit	/	Temp. Coeff. Unit	%/°C
18	Ref. Temp.				20C°, 25C° (68F°, 77F°)
19	Ref. Temp. Unit				°C, °F
20	Cell K				X.X
21	Cell K unit				/cm

Note:

- When printing out dual channel measurement result, conductivity results will be printed after pH results and they will be shown in the same line.
- In continuous mode, “Continuous” will be used for Sample ID.

8.2. Output format

When printed through a printer, the rows of the printing result will be separated and shown as follows:

- The first row: the above item number 1, 2, and 3
- The second row: the above item number 4, 5, 6 and 7
- The third row: the above item number 8, 9, 10, 11 and 12
- The fourth row: the above item number 13, 14 and 15
- The fifth row: the above item number 16, 17, 18 and 19
- The sixth row: the above item number 20 and 21

9. Commands

Commands listed in the following table will be acknowledged by the meter/analyzer.

Table 9-1 Commands

Command Characters	Parameter	Function
SP	/	Print when the measurement is stable
LP	LP 0	Print calibration results
	LP1	Print measurement results
PV	/	Print version number
PSN	/	Print SN number
ON	/	Turn on the meter/analyzer
OFF	/	Turn off the meter/analyzer

10. Maintenance

10.1. Error message

a-AB33PH and a-AB33M1 models:

Table 10-1 Error message

Error Code	Reasons	Solution
Error 3	Measured buffer temperature out of range ($< 0^{\circ}\text{C}$ or $> 50^{\circ}\text{C}$).	Keep the pH buffer temperature within the range for calibration.
Error 4	Offset out of range offset $> 60\text{mV}$ or $< -60\text{mV}$	• Make sure the pH buffer is correct and fresh. • Clean or replace the pH electrode.
Error 5	Slope out of range	• Make sure the buffer is correct and fresh. • Clean or replace the pH electrode.
Error 6	Meter cannot recognize the buffer	• Make sure the buffer is correct and fresh. • Check if the buffer has not been used more than once.

a-AB33EC and a-AB33M1 models:

Table 10-2 Error message

Error Code	Reasons	Solution
Error 3	Measured standard temperature out of range ($< 0^{\circ}\text{C}$ or $> 35^{\circ}\text{C}$).	Keep the standard temperature within the range for calibration
Error 7	Cell K value out of range	• Make sure the Cell K value set in the menu is in compliance with the cell K value of the electrode you use. • The electrode may be broken or expired, so replace the electrode.

Note: If an error happens, the meter will **beep 2 times** to alert.

If the troubleshooting section does not resolve your problem, contact an OHAUS technical support agent.

For assistance in the United States, call toll-free 1-800-526-0659 between 8:00 AM and 5:00 PM Eastern Standard Time. An OHAUS Technical Support Specialist will be available to assist you.

Outside the USA, please visit our website www.ohaus.com to locate the OHAUS office nearest you.

10.2. Meter maintenance

Never unscrew the two halves of the housing!

The meter does not require any maintenance other than occasional wipe with a damp cloth and the replacement of depleted batteries.

The housing is made of acrylonitrile butadiene styrene (ABS). This material does not hold up well to some organic solvents such as toluene, xylene and methyl ethyl ketone (MEK). Any spillage should be immediately wiped off.

10.3. Electrode maintenance

10.3.1 pH electrode maintenance

pH electrode storage

Make sure the electrode is filled with electrolyte solution. Always store the electrode according to the electrode instruction manuals and do not allow it to dry out. It is necessary to store a pH electrode in a solution in order to keep the glass membrane of the pH electrode hydrated. Ideally a storage solution should be used. Never store an electrode in distilled or deionized water.

If the electrode response becomes sluggish or the slope is not acceptable, try the following:

- Soak the electrode in 0.1M HCl for at least 8 hours.
- For fat or oil contaminant, degrease the membrane with cotton wool soaked in either acetone or a soap solution. After electrode treatment, a new calibration should be performed. If the electrode slope is still not acceptable, the electrode might need to be replaced.

10.3.2 Conductivity electrode maintenance

Conductivity electrode storage

Conductivity electrodes can be stored in distilled or deionized water between measurements. For overnight storage or long-term storage, conductivity cells should be thoroughly rinsed and stored dry.

10.4. Cleaning

Warning: Electric Shock Hazard. Disconnect the equipment from the power supply before cleaning.
Make sure that no liquid enters the interior of the instrument.

Attention: Do not use solvents, harsh chemicals, ammonia or

The housing may be cleaned with a cloth dampened with a mild detergent if necessary.

Table 10-3 Cleaning

Pollutant	Cleaning solution	Cleaning time
Water-soluble pollutant	Rinse with deionized water	No limit

Lubricants and oil	Immerse in warm water at about 60°C and washing liquid	10 to 30 minutes
Lime or hydroxide coating	Immerse in 10% acetic acid or 10% hydrochloric acid.	10 to 30 minutes

10.5. Customer service

If you have any questions or need assistance, please contact our technical service experts: Email: pH@Ohaus.com

For other product information, please contact local authorized distributors, local technical sales representatives or contact us using the contact information at the back of this manual.

Visit www.ohaus.com for more products and download product materials, user's guides and manuals, software updates and other application and technical resources.

11. Technical data

11.1. Specification

Equipment Ratings

Indoor use only

Altitude: Up to 2000 m

Operating temperature: 5-40°C

Humidity: Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40°C.

Electrical supply: 12VDC, 0.5A. (For use with certified or approved power supply, which must have a SELV and limited energy output.)

Voltage fluctuation: Mains supply voltage fluctuations up to $\pm 10\%$ of the nominal voltage.

Overvoltage category (installation category): II

Pollution degree: 2

Table 11-1 Specification 1

Benchtop Model		a-AB33PH, a-AB33EC	a-AB33M1
Display Type		6.5" Segment & Dot matrix LCD with backlight	
Multilanguage Interface	User	English, Spanish, French, Portuguese, Chinese, Russian and Turkish	
Measurement Channels		1	2
Measurement Modes	End-point	Auto-stop, Continuous	
Result Log		1000 sets data points with time and date stamp	
Cal. Log		Last calibration	
Keypad		Capacitive touch	
Output		Connect to PC via RS232 and USB, Connect to Printer via RS232.	
pH Electrode Input		BNC	
Conductivity Input		Mini-Din	
Temperature Input		Cinch, NTC 30 kΩ	
Installation Overvoltage		Category II	
Pollution		Degree 2	
Power Supply		Power input: 100-240V~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A	
Net Weight		2.5 kg/5.5lb	
Gross Weight		2.8 kg/7.2lb	
Transportation Size		370 × 270 × 250 mm/14.56 × 10.63 × 9.84 in	

Table 11-2 Specification 2

pH Meter Model		a-AB33PH
Measurement Channels		pH/mV with temperature
PH	Measurement Range	2.00 to 16.00 pH
	Resolution	0.1/0.01 pH
	Resolution Switching	Available
	Accuracy	± 0.01 pH
	Pre-Defined Buffer Groups	3
ORP, RmV	Measurement Range	± 2000.0 mV
	Resolution	1mV
	Accuracy	± 1 mV
	Units	mV, RmV
Temperature	Measurement Range	-5.0 to 110°C
		23°F to 230°F
	Resolution	0.1°C, 0.1°F
	Accuracy	$\pm 0.5^\circ\text{C}$, $\pm 0.5^\circ\text{F}$
	Calibration	None
Calibration	Calibration points	Up to 3 points
	Calibration Sign	Slope/offset & face
	Calibration Mode	Linear

Table 11-3 Specification 3

Conductivity Meter Model		a-AB33EC
Measurement Channels		Conductivity/TDS/Salinity/Resistivity with temperature
Conductivity	Measurement Range	0.01 $\mu\text{S/cm}$ to 19.99 $\mu\text{S/cm}$ 20 $\mu\text{S/cm}$ to 199.9 $\mu\text{S/cm}$ 200 $\mu\text{S/cm}$ to 1999 $\mu\text{S/cm}$ 2.00 mS/cm to 19.99 mS/cm 20.0 mS/cm to 199.9 mS/cm
	Resolution	0.01 $\mu\text{S/cm}$ minimum, auto ranging
	Accuracy	$\pm 0.5\%$ Reading ± 2 LSD
	Ref. Temp.	20 $^{\circ}\text{C}$, 25 $^{\circ}\text{C}$
	Cell K	0.01 to 10.00 cm^{-1}
	Temperature Compensation	Linear (0 to 10.0%/ $^{\circ}\text{C}$), off
TDS	Measurement Range	0.1 mg/L to 200 g/L
	Resolution	0.01 mg/L minimum, auto ranging
	Accuracy	$\pm 0.5\%$ Reading ± 2 LSD
	TDS Factor Range	Linear, 0.04 to 10.00, default 0.5
Resistivity	Measurement Range	1 to 100 $\text{M}\Omega\text{-cm}$
	Resolution	0.01 $\Omega\text{-cm}$, auto ranging
	Accuracy	$\pm 0.5\%$ Reading ± 2 LSD
Practical SALT	Measurement Range	0 to 100 psu
	Resolution	0.01 psu minimum, auto ranging
	Accuracy	$\pm 0.5\%$ Reading ± 2 LSD
Temperature	Measurement Range	-5.0 to 110.0 $^{\circ}\text{C}$, 32.0 $^{\circ}\text{F}$ to 212.0 $^{\circ}\text{F}$
	Resolution	0.1 $^{\circ}\text{C}$, 0.1 $^{\circ}\text{F}$
	Accuracy	$\pm 0.3^{\circ}\text{C}$, $\pm 0.3^{\circ}\text{F}$
	Calibration	None
Calibration	Calibration points	Calibration of one-point Cell K; 6 available buffers (10 $\mu\text{S/cm}$, 84 $\mu\text{S/cm}$, 146.5 $\mu\text{S/cm}$, 500 $\mu\text{S/cm}$, 1413 $\mu\text{S/cm}$, 12.88 mS/cm).
	Calibration Sign	Cell constant & face
	Calibration Mode	Linear

Table 11-4 Specification 4

pH and conductivity analyzer model		a-AB33M1
Measurement Channels		Change of pH/mV/ORP/conductivity/TDS/SALT/resistivity with temperature
pH	Measurement Range	-2.00 to 20.00 pH
	Resolution	0.01 pH
	Resolution Switching	Available
	Accuracy	± 0.01 pH
	Pre-Defined Buffer Groups	5
ORP. Redox	Measurement Range	± 2000.0 mV
	Resolution	0.1mV
	Accuracy	± 0.5 mV
	Units	mV, RmV
Conductivity	Measurement Range	0.01 μ S/cm to 19.99 μ S/cm 20 μ S/cm to 199.9 μ S/cm 200 μ S/cm to 1999 μ S/cm 2.00mS/cm to 19.99mS/cm 20.0mS/cm to 500.0mS/cm
	Resolution	0.01 μ S/cm minimum, auto ranging
	Accuracy	± 0.5 % Reading ± 2 LSD
	Ref. Temp.	20°C, 25°C
	Cell K	0.001 to 10.000 cm ⁻¹
	Temperature Compensation	Linear (0 to 10.0%/°C), off
TDS	Measurement Range	0.1 mg/L to 199.9 g/L
	Resolution	0.01 mg/L minimum, auto ranging
	Accuracy	± 0.5 % Reading ± 2 LSD
	TDS Factor Range	Linear, 0.01 to 10.00, default 0.5
Resistivity	Measurement Range	2 to 100 M Ω -cm
	Resolution	0.01 Ω -cm, auto ranging
	Accuracy	± 0.5 % Reading ± 2 LSD
Practical SALT	Measurement Range	0 to 100 psu
	Resolution	0.01 psu minimum, auto ranging
	Accuracy	± 0.5 % Reading ± 2 LSD
Temperature	Measurement Range	-5.0 to 110°C, 23°F to 230°F
	Resolution	0.1°C, 0.1°F

	Accuracy	$\pm 0.3^{\circ}\text{C}$, $\pm 0.3^{\circ}\text{F}$
	Calibration	None
Calibration	Calibration points	Up to 5 points for pH; 1 point cell constant calibration; 6 buffers available (10, 84, 146.5, 500, 1413 $\mu\text{S}/\text{cm}$, 12.88 mS/cm)
	Calibration Sign	Slope/offset & Cell constant & face
	Calibration Mode	Linear

11.2. Size


Table 11-5 Size

12. Appendixes

12.1. pH buffer group

The meter/analyzer automatically corrects for the temperature dependence of the buffer group of the pH value given in the following tables.

Table 12-1 Buffer group: USA standard

Temperature °C	pH1.68	pH4.01	pH7.00	pH10.00
5	1.67	4.01	7.09	10.25
10	1.67	4.00	7.06	10.18
15	1.67	4.00	7.04	10.12
20	1.68	4.00	7.02	10.06
25	1.68	4.01	7.00	10.01
30	1.68	4.01	6.99	9.97
35	1.69	4.02	6.98	9.93

Table 12-1 Buffer group: JJG standard

Temperature °C	pH4	pH6.86	pH9.18	pH12.45
5	3.99	6.95	9.40	12.73
10	4.00	6.93	9.33	12.67
15	4.00	6.89	9.27	12.59
20	4.00	6.88	9.23	12.52
25	4.00	6.86	9.18	12.45
30	4.01	6.85	9.14	12.36
35	4.02	6.84	9.10	12.25

Table 12-2 Buffer group: European standard

Temperature °C	pH2.00	pH4.01	pH7.00	pH9.21	pH11.00
5	2.02	4.01	7.09	9.45	11.72
10	2.01	4.00	7.06	9.38	11.54
15	2.00	4.00	7.04	9.32	11.36
20	2.00	4.00	7.02	9.26	11.18
25	2.00	4.01	7.00	9.21	11.00
30	1.99	4.01	6.99	9.16	10.82
35	1.99	4.02	6.98	9.11	10.64
40	1.98	4.03	6.97	9.06	10.46
45	1.98	4.04	6.97	9.03	10.28
50	1.98	4.06	6.97	8.99	10.10

Table 12-3 Buffer group: MERCK (Ref. Temp. 20°C)

Temperature °C	pH2.00	pH4.00	pH7.00	pH9.00	pH12.00
5	2.01	4.04	7.07	9.16	12.41
10	2.01	4.02	7.05	9.11	12.26
15	2.00	4.01	7.02	9.05	12.10
20	2.00	4.00	7.00	9.00	12.00
25	2.00	4.01	6.98	8.95	11.88
30	2.00	4.01	6.98	8.91	11.72
35	2.00	4.01	6.96	8.88	11.67
40	2.00	4.01	6.95	8.85	11.54
45	2.00	4.01	6.95	8.82	11.44
50	2.00	4.00	6.95	8.79	11.33

Table 12-4 Buffer group: DIN (19266) (Ref. Temp. 25°C)

Temperature °C	pH1.09	pH4.65	pH6.79	pH9.23	pH12.75
5	1.668	4.004	6.950	9.392	13.207
10	1.670	4.001	6.922	9.331	13.003
15	1.672	4.001	6.900	9.277	12.810
20	1.676	4.003	6.880	9.228	12.627
25	1.680	4.008	6.865	9.184	12.454
30	1.685	4.015	6.853	9.144	12.289
35	1.691	4.026	6.845	9.110	12.133
40	1.697	4.036	6.837	9.076	11.984
45	1.704	4.049	6.834	9.046	11.841
50	1.712	4.064	6.833	9.018	11.705

12.2. Conductivity standards

Table 12-5 Conductivity standard

Temperature °C	10 µS/cm	84 µS/cm	146.5 µS/cm	500µS/cm	1413 µS/cm	12.88 mS/cm
5	6.1 µS/cm	53 µS/cm	93 µS/cm	315 µS/cm	896 µS/cm	8.22 mS/cm
10	7.1 µS/cm	60 µS/cm	106 µS/cm	360 µS/cm	1020 µS/cm	9.33 mS/cm
15	8.0 µS/cm	68 µS/cm	119 µS/cm	403 µS/cm	1147 µS/cm	10.48 mS/cm
20	9.0 µS/cm	76 µS/cm	133 µS/cm	452 µS/cm	1278 µS/cm	11.67 mS/cm

25	10.0 μS/cm	84 μS/cm	146.5 μS/cm	500 μS/cm	1413 μS/cm	12.88 mS/cm
30	11.0 μS/cm	92 μS/cm	161 μS/cm	549 μS/cm	1552 μS/cm	14.12 mS/cm
35	12.1 μS/cm	101 μS/cm	176 μS/cm	603 μS/cm	1667 μS/cm	15.39 mS/cm

12.3. Example of Temp. Coeff. (α value)

Table 12-6 Example of Temp. Coeff.

Substance at 25°C	Concentration [%]	Temp. Coeff. α [%/°C]
Hydrogen chloride	10	1.56
Potassium chloride	10	1.88
Acetic acid	10	1.69
Sodium chloride	10	2.14
Sulfuric acid	10	1.28
Hydrofluoric acid	1.5	7.2

α -coefficient of conductivity standard used to calculate the Ref. Temp. of 25°C.

Table 12-7 α coefficient of conductivity standard

Standard	Measurement Temperature: 15°C	Measurement Temperature: 20°C	Measurement Temperature: 30°C	Measurement Temperature: 35°C
84μS/cm	1.95	1.95	1.95	2.01
1413μS/cm	1.94	1.94	1.94	1.99
12.88mS/cm	1.90	1.89	1.91	1.95

13. Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Marking	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC), 2014/35/EU (LVD). The EU Declaration of Conformity is available online at www.ohaus.com/ce .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.ohaus.com/wEEE .
	EN 61326-1

ISED Canada Compliance Statement:

CAN ICES-003 (A) / NMB-003 (A)

International organization for standardization (ISO) 9001 registration

The management system governing the production of this product is ISO 9001 certified.

FCC Supplier Declaration of Conformity

Unintentional Radiator per 47CFR Part B
Trade Name: OHAUS CORPORATION
Model or Family identification: Aquasercher a-AB

Issuing Party that Assembled the Product:

Ohaus Instruments (Changzhou) Co., Ltd.
Building C, No. 6 Zhengqiang Road, Xuejia Town, Xinbei District, Changzhou
Jiangsu 213022
China
Phone: +86 519 85287270

Responsible Party – U.S. Contact Information:

Ohaus Corporation
8 Campus Drive, Suite 105
Parsippany, NJ 07054
United States
Phone: +1 973 377 9000
Web: www.ohaus.com

FCC Compliance Statement:

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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OHAUS® and Aquasearcher are either registered trademarks or trademarks of OHAUS Corporation in the United States and/or other countries.
P/N 31012761 A © 2024 OHAUS Corporation, all rights reserved.
OHAUS Corporation